 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Loading
	Test Case Number	LIMS-SM-012
	Actual Result	Pass

1. Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001





The screenshot shows the LIMS login interface. At the top center, the word "LIMS" is displayed in yellow. Below it is the PharmaSoftSol logo. The login form contains three input fields: "User Name" with the value "LIMS_APPROVER", "Password" with masked characters "*****", and "Module" with a dropdown menu set to "STABILITY MANAGEMENT". At the bottom of the form are two buttons: "Login" and "Reset".

1. Click on Loading

SSNo : PSS-LIMS-SS-001

PHARMA SOFT
STABILITY MANAGEMENT



LIMS APPROVER

Loading Approval

Activity :

- . Select The Record and Click On View.
- . Enter Comments and Click on Approve / Revert.
- . Enter E-Signature Password and Click on Submit.

S.No.	Protocol Number	Product Name	Generic Name	Specification Name	Loading By/On	Details
1	PS0001	Aspirin	Aspirin and codeine	SS Aspirin Spec	LIMS USER 09/06/2022 20:18:09	View
2	PS0002	Aspirin	Aspirin and codeine	SS Aspirin Spec	LIMS USER 09/06/2022 20:53:37	View
3	PS0024	Aspirin	Aspirin and codeine	SS Aspirin Spec	LIMS USER 15/07/2022 11:00:50	View

<< < 1 > >>
Go to page 1
Row count: 10
Showing 1-3 of 3

LIMS 2.0
All © rights Reserved PharmaSoftSol

2. Select a Record And Click on View

SSNo : PSS-LIMS-SS-002

Test Case Name	Approve Loading
Test Case Number	LIMS-SM-012
Actual Result	Pass

Loading Approval
Loading Details ✕

Activity :
 . Select The Record at
 . Enter Comments and
 . Enter E-Signature Pa

Loading Details

S.No.	Protocol
1	PS0001
2	PS0002
3	PS0024

<< < 1 > >> Go to p


Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Box	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	15/07/2022	<button>View</button>

Comments*:

Approve
Revert

Details
 20:18:09 View
 20:53:37 View
 11:00:50 View
 Showing 1-3 of 3

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Loading
	Test Case Number	LIMS-SM-012
	Actual Result	Pass

3. Enter Comments

SSNo : PSS-LIMS-SS-003

PHARMA SOFT

STABILITY MANAGEMENT

LIMS APPROVER

Loading Approval

Activity :
. Select The Record at
. Enter Comments and
. Enter E-Signature Pa

Loading Details

S.No.	Protocol
1	PS0001
2	PS0002
3	PS0024

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Box	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	15/07/2022	View

Comments * :

Approved

Approve

Revert

Details

20:18:09 [View](#)
20:53:37 [View](#)
11:00:50 [View](#)
Showing 1-3 of 3

LIMS 2.0
All © rights Reserved PharmaSoftSol

4. Click on Approve

SSNo : PSS-LIMS-SS-004

PHARMA SOFT

STABILITY MANAGEMENT

LIMS APPROVER

Loading Approval

Activity :
. Select The Record at
. Enter Comments and
. Enter E-Signature Pa

Loading Details

S.No.	Protocol
1	PS0001
2	PS0002
3	PS0024

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Box	44	40°C ± 2°C / 75%RH ± 5%RH	UP R					Gm			Gm	15/07/2022	View

Comments * :

Appro

Validate E-Signature

E-Signature Password* :


Submit

Cancel

Details

20:18:09 [View](#)
20:53:37 [View](#)
11:00:50 [View](#)
Showing 1-3 of 3

LIMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Loading
	Test Case Number	LIMS-SM-012
	Actual Result	Pass

5. Enter E-Signature Password

SSNo : PSS-LIMS-SS-005

PHARMA SOFT **STABILITY MANAGEMENT** ☰ 👤 LIMS APPROVER

Loading Approval

Activity :
 . Select The Record at
 . Enter Comments and
 . Enter E-Signature Pa

Loading Details

S.No.	Protocol
1	PS0001
2	PS0002
3	PS0024

<< < 1 > >> Go to p

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Box	44	40°C ± 2°C / 75%RH ± 5%RH	UP R								Gm	15/07/2022	View

Comments * : Appro

Validate E-Signature

E-Signature Password* :

[Submit](#) [Cancel](#)

Details

20:18:09	View
20:53:37	View
11:00:50	View

Showing 1-3 of 3

6. Click on Submit

SSNo : PSS-LIMS-SS-006

PHARMA SOFT **STABILITY MANAGEMENT** ☰ 👤 LIMS APPROVER

Loading Approval

Activity :
 . Select The Record at
 . Enter Comments and
 . Enter E-Signature Pa

Loading Details

S.No.	Protocol
1	PS0001
2	PS0002
3	PS0024

<< < 1 > >> Go to p

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Box	44	40°C ± 2°C / 75%RH ± 5%RH	UP R								Gm	15/07/2022	View

Comments * : Appro

Validate E-Signature


E-Signature Password* :

[Saving...](#) [Submit](#) [Cancel](#)

Details

20:18:09	View
20:53:37	View
11:00:50	View

Showing 1-3 of 3

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Loading
	Test Case Number	LIMS-SM-012
	Actual Result	Pass

7. Click on OK button

SSNo : PSS-LIMS-SS-007

PHARMA SOFT
STABILITY MANAGEMENT
☰
 LIMS APPROVER ▾

Loading Approval

Activity :

- . Select The Record and
- . Enter Comments and
- . Enter E-Signature Pa

Loading Details

S.No.	Protocol
1	PS0001
2	PS0002
3	PS0024

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/F Details
Box	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	15/07/2022	View

Comments : Approved

Success

Submitted successfully.

OK

Approve

Revert

Details

20:18:09 [View](#)
20:53:37 [View](#)
11:00:50 [View](#)

Showing 1-3 of 3

LIMS 2.0

All © rights Reserved PharmaSoftSol

8. Click on User Name


SSNo : PSS-LIMS-SS-008

PHARMA SOFT
STABILITY MANAGEMENT
☰
 LIMS APPROVER ▾

Loading Approval

Activity :

- . Select The Record and Click On View.
- . Enter Comments and Click on Approve / Revert.
- . Enter E-Signature Password and Click on Submit.



Loading Details


S.No.	Protocol Number	Product Name	Generic Name	Specification Name	Loading By/On	Logout	Details
1	PS0001	Aspirin	Aspirin and codeine	SS Aspirin Spec	LIMS USER 09/06/2022 20:18:09		View
2	PS0002	Aspirin	Aspirin and codeine	SS Aspirin Spec	LIMS USER 09/06/2022 20:53:37		View

Go to page: 1 Row count: 10

Showing 1-2 of 2

LIMS 2.0

All © rights Reserved PharmaSoftSol

	Test Case Name	Approve Loading
	Test Case Number	LIMS-SM-012
	Actual Result	Pass

9. Click on Logout

SSNo : PSS-LIMS-SS-009

