	<b>Test Case Name</b>	Approve Loading
	<b>Test Case Number</b>	LIMS-SM-012
	<b>Actual Result</b>	Pass

1. Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001



1. Click on Loading


SSNo : PSS-LIMS-SS-001



2. Select a Record And Click on View

SSNo : PSS-LIMS-SS-002

LIMS 2.0

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Approve Loading
	<b>Test Case Number</b>	LIMS-SM-012
	<b>Actual Result</b>	Pass

### 3. Enter Comments

SSNo : PSS-LIMS-SS-003

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver

Loading Approval

Activity :  
. Select The Record at  
. Enter Comments and  
. Enter E-Signature Pa

Loading Details  
S.No. Protocol  
1 PC0012  
<< < 1 > >> Go to p

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	2	3	Gm	3	2	Gm	31/08/2022	<a href="#">View</a>

Comments\* :

Approved

Approve
Revert

Details  
6:51:35 [View](#)  
Showing 1-1 of 1

LIMS 2.0
All © rights Reserved PharmaSoftSol

### 4. Click on Approve

SSNo : PSS-LIMS-SS-004

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver

Loading Approval

Activity :  
. Select The Record at  
. Enter Comments and  
. Enter E-Signature Pa

Loading Details  
S.No. Protocol  
1 PC0012  
<< < 1 > >> Go to p

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP R								Gm	31/08/2022	<a href="#">View</a>

Comments\* :

Appro


Validate E-Signature

E-Signature Password\* :

Submit
Cancel

Details  
6:51:35 [View](#)  
Showing 1-1 of 1

LIMS 2.0
All © rights Reserved PharmaSoftSol

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Approve Loading
	<b>Test Case Number</b>	LIMS-SM-012
	<b>Actual Result</b>	Pass

## 5. Enter E-Signature Password

SSNo : PSS-LIMS-SS-005

**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS Approver

Loading Approval

Activity :  
 . Select The Record at  
 . Enter Comments and  
 . Enter E-Signature Pa

**Loading Details**

S.No. Protocol  
 1 PC0012

<< < 1 > >> Go to Page

**Loading Details**

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Reference Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP R					Gm				31/08/2022	<a href="#">View</a>

Comments\* : Appro

**Validate E-Signature**

E-Signature Password\* :

[Submit](#) [Cancel](#)

LIMS 2.0 All © rights Reserved PharmaSoftSol

## 6. Click on Submit

SSNo : PSS-LIMS-SS-006

**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS Approver

Loading Approval

Activity :  
 . Select The Record at  
 . Enter Comments and  
 . Enter E-Signature Pa

**Loading Details**

S.No. Protocol  
 1 PC0012

<< < 1 > >> Go to Page


**Loading Details**

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Reference Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	2	3	Gm	3	2	Gm	31/08/2022	<a href="#">View</a>

Comments\* : Approved

[Approve](#) [Revert](#)

LIMS 2.0 All © rights Reserved PharmaSoftSol


	<b>Test Case Name</b>	Approve Loading
	<b>Test Case Number</b>	LIMS-SM-012
	<b>Actual Result</b>	Pass

7.Click on OK button

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT**

**STABILITY MANAGEMENT**

 LIMS Approver

Loading Approval

Activity :  
. Select The Record and  
. Enter Comments and  
. Enter E-Signature Pa

Loading Details

S.No. Protocol  
1 PC0012

<< < 1 > >> Go to

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/F Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	2	3	Gm	3	2	Gm	31/08/2022	<a href="#">View</a>

Comments : Approved

Success

Submitted successfully.

OK

Approve Revert

Details

6:51:35 View

Showing 1-1 of 1

LIMS 2.0

All © rights Reserved PharmaSoftSol

8.Click on User Name

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT**

**STABILITY MANAGEMENT**

 LIMS Approver

Loading Approval


Activity :  
. Select The Record and Click On View.  
. Enter Comments and Click on Approve / Revert.  
. Enter E-Signature Password and Click on Submit.

Loading Details

S.No. Protocol Number Product Name Generic Name Specification Name Loading By/On Logout Details


Row count: 10

No data available!



LIMS 2.0

All © rights Reserved PharmaSoftSol

	<b>Test Case Name</b>	Approve Loading
	<b>Test Case Number</b>	LIMS-SM-012
	<b>Actual Result</b>	Pass

9. Click on Logout

SSNo : PSS-LIMS-SS-009

