
	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

1.Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001



The screenshot shows the LIMS login interface. It features the PharmaSoftSol logo at the top left. The login form includes fields for 'User Name' (containing 'LIMS_USER'), 'Password' (masked with dots), and 'Module' (a dropdown menu set to 'STABILITY MANAGEMENT'). Below these fields are 'Login' and 'Reset' buttons. The background is a blue-to-purple gradient.


2.Click on Sample Acknowledgement

SSNo : PSS-LIMS-SS-002



The screenshot displays the 'Sample Acknowledgement' section of the LIMS application. At the top, there's a header with 'PHARMA SOFT' and 'STABILITY MANAGEMENT'. Below this, the 'Sample Acknowledgement' title is followed by an 'Action' button and a 'Sample Acknowledgment Details' link. Instructions for selecting and reverting sample acknowledgments are provided. Two radio buttons are present: 'Sample Acknowledgment' (unselected) and 'Reverted Sample Acknowledgments by Approver' (selected). Below the radio buttons is a table titled 'List of Reverted Sample Acknowledgements by Approver'. The table has columns for S.No, Product Name, Specification Number, AR Number, Protocol Number, Storage Conditions, Orientation, Stability Conditions, Packing Type, Batch Number, Start Date, RevertedBy/On, and Comments. A single record is shown for product 'PARA6S' with a 'view' button in the comments column. At the bottom, there are pagination controls showing 'Go to page: 1' and 'Row count: 10', along with 'Showing 1-1 of 1'.

S.No	Product Name	Specification Number	AR Number	Protocol Number	Storage Conditions	Orientation	Stability Conditions	Packing Type	Batch Number	Start Date	RevertedBy/On	Comments
1	PARA6S	SPSSPARA6S	ASS220051	PC0007	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	Blister	Product Batch No 001	17/08/2022	LIMS User 17/08/2022	view

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

3. Click on Sample Acknowledgment

SSNo : PSS-LIMS-SS-003

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action Sample Acknowledgment Details

. Select Sample Acknowledgment, Click On Search, Select The Record and Click on Next.
. If Reverted Sample Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.

☒ Sample Acknowledgment
☐ Reverted Sample Acknowledgments by Approver

From Date : To Date:
Product : **Browse** Specification:

Search **Clear**

Samples Loaded														
S.No	Product	Specification Number	AR Number	Protocol Number	Storage Condition	Stability Condition	Orientation	Packing Type	Batch Number	Start Date	Total Qty Remained	Qty Remained UOM	Interval	Is deviated
No data available!														

Go to page: Row count:

IMS 2.0
All © rights Reserved PharmaSoftSol

4. Select a Record

SSNo : PSS-LIMS-SS-004

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action Sample Acknowledgment Details

. Select Sample Acknowledgment, Click On Search, Select The Record and Click on Next.
. If Reverted Sample Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.


☒ Sample Acknowledgment
☐ Reverted Sample Acknowledgments by Approver

From Date : To Date:
Product : **Browse** Specification:

Search **Clear**

Samples Loaded														
S.No	Product	Specification Number	AR Number	Protocol Number	Storage Condition	Stability Condition	Orientation	Packing Type	Batch Number	Start Date	Total Qty Remained	Qty Remained UOM	Interval	Is deviated
1	Epinephrine	SS SN 001	ASS220001	PC0002	40°C ± 2°C / 75%RH ± 5%RH 30°C ± 2°C	ACCELERATED	SIDE BY SIDE	Blister	BN 001	12/08/2022	2	mg / mL	1 DAY	No

IMS 2.0
All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

5. Click on Next

SSNo : PSS-LIMS-SS-005

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

[Action](#)
[Sample Acknowledgment Details](#)

On Search, Select The Record and Click on Next.
by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.

☐ Reverted Sample Acknowledgments by Approver

To Date:

Specification:

Number	Protocol Number	Storage Condition	Stability Condition	Orientation	Packing Type	Batch Number	Start Date	Total Qty Remained	Qty Remained UOM	Interval	Is deviated
1001	PC0002	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	SIDE BY SIDE	Blister	BN 001	12/08/2022	2	mg / mL	1 DAY	No
		30°C ± 2°C									

IMS 2.0
All © rights Reserved PharmaSoftSol

6. Enter Quantity

SSNo : PSS-LIMS-SS-006

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

[Action](#)
[Sample Acknowledgment Details](#)

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.


Quantity Withdraw :
UOM :
Approval From:

Used For:

Maximum allowed characters should be 100

Select an Action*:
☐ Acknowledge
☐ Terminate

IMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

7.Select UOM

SSNo : PSS-LIMS-SS-007

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action
Sample Acknowledgment Details

- Enter Quantity withdraw and select UOM from drop down.click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :
UOM :
Approval From:

Used For:

Select an Action*: ☐ Acknowledge ☐ Terminate

IMS 2.0
All © rights Reserved PharmaSoftSol

8.Enter Used For

SSNo : PSS-LIMS-SS-008

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action
Sample Acknowledgment Details


- Enter Quantity withdraw and select UOM from drop down.click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :
UOM :
Approval From:

Used For:


Select an Action*: ☐ Acknowledge ☐ Terminate

IMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

11.Select a Record

SSNo : PSS-LIMS-SS-011




The screenshot shows the PHARMA SOFT STABILITY MANAGEMENT interface. The 'Locations' section on the left has a tree view with 'organization' and 'Plant-1' (containing QA, QC, and EN). The 'Users' section on the right has a table with columns 'First Name' and 'Last Name'. The table contains 5 rows: LIMS User, LIMS Approver, LIMS Reviewer, User 01, and User 02. The 'LIMS Approver' row is selected. Below the table are pagination controls showing 'Showing 1-5 of 12' and buttons for 'Select' and 'Cancel'.

12.Click on Select

SSNo : PSS-LIMS-SS-012



The screenshot shows the PHARMA SOFT STABILITY MANAGEMENT interface with the 'Sample Acknowledgment' form. The form includes fields for 'Quantity Withdraw' (set to 2), 'UOM' (set to Gm), and 'Approval From' (set to LIMS Approver). The 'Used For' field contains 'Sampling'. Below these fields are radio buttons for 'Acknowledge' and 'Terminate'. The 'LIMS Approver' button is highlighted with a 'Select' label.

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

13. Click on Acknowledge

SSNo : PSS-LIMS-SS-013

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action
Sample Acknowledgment Details

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :
UOM :
Approval From:

Used For:

Select an Action*: ☒ Acknowledge ☐ Terminate

LIMS 2.0
All © rights Reserved PharmaSoftSol

14. Click on Finish

SSNo : PSS-LIMS-SS-014

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action
Sample Acknowledgment Details

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.


Quantity Withdraw :
UOM :
Approval From:

Used For:

Select an Action*: ☒ Acknowledge ☐ Terminate


Validate E-Signature
E-Signature Password*:

LIMS 2.0
All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

15. Enter E-Signature Password

SSNo : PSS-LIMS-SS-015



PHARMA SOFT **STABILITY MANAGEMENT** IMS User

Sample Acknowledgment

Action **Sample Acknowledgment Details**

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments, or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw : 2

Used For: Sampling

Select an Action*: ☒ Acknowledgment ☐ Terminate

IMS Approver Select

Validate E-Signature


E-Signature Password*:

Submit Cancel

IMS 2.0 All © rights Reserved PharmaSoftSol

16. Click on Submit

SSNo : PSS-LIMS-SS-016



PHARMA SOFT **STABILITY MANAGEMENT** IMS User

Sample Acknowledgment

Action **Sample Acknowledgment Details**

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments, or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw : 2

Used For: Sampling

Select an Action*: ☒ Acknowledgment ☐ Terminate


IMS Approver Select

Validate E-Signature

E-Signature Password*:

Saving... Submit Cancel

IMS 2.0 All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

17. Click on OK button

SSNo : PSS-LIMS-SS-017

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action **Sample Acknowledgment Details**

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :
UOM :

Used For:

Select an Action*: ☒ Acknowledge

Success

Submitted successfully.

OK

LIMS Approver
Select

LIMS 2.0
All © rights Reserved PharmaSoftSol

18. Click on User Name

SSNo : PSS-LIMS-SS-018

PHARMA SOFT
STABILITY MANAGEMENT
LIMS User

Sample Acknowledgment

Action **Sample Acknowledgment Details**

- . Select Sample Acknowledgment, Click On Search, Select The Record and Click on Next.
- . If Reverted Sample Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.


☐ Sample Acknowledgment
☒ Reverted Sample Acknowledgments by Approver

S.No	Product Name	Specification Number	AR Number	Protocol Number	Storage Conditions	Orientation	Stability Conditions	Packing Type	Batch Number	Start Date	RevertedBy/On	Comments
1	PARA65	SPSSPARA65	ASS220051	PC0007	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	Blister	Product Batch No 001	17/08/2022	LIMS User 17/08/2022	View

Go to page:
Row count:
Showing 1-1 of 1


Logout

LIMS 2.0
All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

19. Click on Logout

SSNo : PSS-LIMS-SS-019

