	<b>Test Case Name</b>	Reject Protocol
	<b>Test Case Number</b>	PSS-LIMS-005
	<b>Actual Result</b>	Pass

1. Enter the username, password, Select Module and click on login button


SSNo : PSS-LIMS-SS-001



The image shows the LIMS login interface. It features the PharmaSoftSol logo at the top left. The main heading is 'LIMS'. Below it, there are three input fields: 'User Name' with the value 'LIMS\_APPROVER', 'Password' with masked characters '\*\*\*\*\*', and 'Module' with a dropdown menu set to 'STABILITY MANAGEMENT'. At the bottom of the form are two buttons: 'Login' and 'Reset'.

2. Click on Protocol

SSNo : PSS-LIMS-SS-002




The image shows the 'Protocol Approval' page in the LIMS system. The header includes the PharmaSoftSol logo, the text 'STABILITY MANAGEMENT', and a user profile icon labeled 'LIMS Approver'. The main content area is titled 'Protocol Approval' and contains instructions for the activity: 'Select The Record and Click On View & Approval.', 'Enter Comments and Click on Approve / Revert.', and 'Enter E-Signature Password and Click on Submit.'. Below the instructions is a table titled 'List of Protocols for Approval'.

S.No.	Protocol Number	Product Name	Specification Name	Generic Name	Registered By/On	Request From	Action
1	PC0013	PARA65	SPPARACHEMICAL	Paracetamol 650 mg	LIMS User 31/08/2022 17:25:03	NEW	<a href="#">View &amp; Approval</a>

Below the table, there are navigation controls: '<<' '<' '1' '>' '>>' 'Go to page: 1' 'Row count: 10' and 'Showing 1-1 of 1'.

The footer of the page contains 'LIMS 2.0' and 'All © rights Reserved PharmaSoftSol'.

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Protocol
	<b>Test Case Number</b>	PSS-LIMS-005
	<b>Actual Result</b>	Pass

### 3. Select a Record And Click on View & Approval

SSNo : PSS-LIMS-SS-003

**PHARMA SOFT**
**STABILITY MANAGEMENT**
88
LIMS Approver

rotocol Stability Plan

Procedure :

↶

↷

Formats

B

I

≡

≡

≡

≡

≡

≡

≡

≡

rotocol Procedure

Location:

Plant-1

Specification Number:

SPPARACH001

Generic Name:

Paracetamol 650 mg

Test Planner Version:

0

Product Name:

PARA65

Specification Name:

SPPARACHEMICAL

STP Number:

STPCHPARA001

Customer Name:

BA

### 4. Enter Comments

SSNo : PSS-LIMS-SS-004

**PHARMA SOFT**
**STABILITY MANAGEMENT**
88
LIMS Approver


Packing type	No	Batch size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Type
Blister	44	100	08/31/2022	08/31/2027	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH

Comments\*:

Rejected

Approve

Revert

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Protocol
	<b>Test Case Number</b>	PSS-LIMS-005
	<b>Actual Result</b>	Pass

5. Click on Reject

SSNo : PSS-LIMS-SS-005

**PHARMA SOFT**

**STABILITY MANAGEMENT**

LIMS Approver

Packing type	No	Batch size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Type
Blister	44	100	08/31/2022	08/31/2027	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH

Comments\*:

Rejected

Validate E-Signature

E-Signature Password\*:

Submit
Cancel

Approve
Revert

6. Enter E-Signature Password

SSNo : PSS-LIMS-SS-006

**PHARMA SOFT**

**STABILITY MANAGEMENT**

LIMS Approver

Packing type	No	Batch size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Type
Blister	44	100	08/31/2022	08/31/2027	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH

Comments\*:


Rejected

Validate E-Signature

E-Signature Password\*:

Submit
Cancel

Approve
Revert

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Protocol
	<b>Test Case Number</b>	PSS-LIMS-005
	<b>Actual Result</b>	Pass

7.Click on Submit

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
👤 LIMS Approver ▾

Packing type	No	Batch size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Type
Blister	44	100	08/31/2022	08/31/2027	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH

Comments\*:

Rejected

Approve
Revert

8.Click on OK button

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
👤 LIMS Approver ▾


Packing Type	Product Batch No	Batch Size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Interval Type	Int
Blister	44	100	08/31/2022	08/31/2027	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	

Comments\*:

Rejected

Success
Submitted successfully.
OK


Approve
Revert

	<b>Test Case Name</b>	Reject Protocol
	<b>Test Case Number</b>	PSS-LIMS-005
	<b>Actual Result</b>	Pass

9.Click on User Name


SSNo : PSS-LIMS-SS-009

**PHARMA SOFT SOL**
**STABILITY MANAGEMENT**

 LIMS Approver

### Protocol Approval

**Activity :**  
 . Select The Record and Click On View & Approval.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.



#### List of Protocols for Approval

S.No.	Protocol Number	Product Name	Specification Name	Generic Name	Registered By/On	Request From	Action
No data available!							

Row count: 10

LIMS 2.0
 All © rights Reserved PharmaSoftSol

10.Click on Logout

SSNo : PSS-LIMS-SS-010

LIMS



User Name

Password

Module

Select

Login

Reset