	Test Case Name	Approve Sample Acknowledgement
	Test Case Number	LIMS-SM-014
	Actual Result	Pass

1. Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001





The login form is titled "LIMS" and features the PharmaSoftSol logo. It contains three input fields: "User Name" with the value "LIMS_APPROVER", "Password" with masked characters "*****", and "Module" with a dropdown menu set to "STABILITY MANAGEMENT". Below the fields are "Login" and "Reset" buttons.

2. Click on Sample Acknowledgement

SSNo : PSS-LIMS-SS-002

PHARMA SOFT
STABILITY MANAGEMENT



LIMS APPROVER

Sample Acknowledgment Approval


Activity :

- . Select The Record and Click On view.
- . Enter Comments and Click on Approve / Revert.
- . Enter E-Signature Password and Click on Submit.

S.No	Protocol Number	Product Name	Specification Number	AR Number	Storage Condition	Stability Condition	Packing type	Batch Number	Start Date	Done By/On	Action
1	PS0014	Aspirin	SNSS 001	ASS220038	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	Product Batch No	21/06/2022	LIMS USER /21/06/2022 14:35:27	View
2	PS0014	Aspirin	SNSS 001	ASS220041	30°C ± 2°C / 65%RH ± 5%RH	INTERMEDIATE	Box	Product Batch No	21/06/2022	LIMS USER /21/06/2022 14:40:00	View
3	PS0023	Aspirin	SNSS 001	ASS220068	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	44	15/07/2022	LIMS USER /15/07/2022 10:48:43	View

Go to page: 1 Row count: 10
Showing 1-3 of 3

LIMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Sample Acknowledgement
	Test Case Number	LIMS-SM-014
	Actual Result	Pass

3. Select a Record And Click on View

SSNo : PSS-LIMS-SS-003

PHARMA SOFT
STABILITY MANAGEMENT
LIMS APPROVER

Sample Acknowledgement
 Activity :
 . Select The
 . Enter Com
 . Enter E-Sig
List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number:
 Product Name:

Specification Number:
 AR Number:

Batch Number:
 Packing Type:

Storage Conditions:
 Interval:

Qty Withdraw: UOM:
 Total Qty Remained: UOM:

Start Date:

Used For:

Comments*:

Action

 Showing 1-3 of 3

LIMS 2.0

4. Enter Comments

SSNo : PSS-LIMS-SS-004

PHARMA SOFT
STABILITY MANAGEMENT
LIMS APPROVER

Sample Acknowledgement
 Activity :
 . Select The
 . Enter Com
 . Enter E-Sig
List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number:
 Product Name:

Specification Number:
 AR Number:

Batch Number:
 Packing Type:

Storage Conditions:
 Interval:

Qty Withdraw: UOM:
 Total Qty Remained: UOM:

Start Date:


Used For:

Comments*:

Action


 Showing 1-3 of 3

LIMS 2.0

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Sample Acknowledgement
	Test Case Number	LIMS-SM-014
	Actual Result	Pass

5. Click on Approve

SSNo : PSS-LIMS-SS-005



PHARMA SOFT **STABILITY MANAGEMENT** LIMS APPROVER

Sample Acknowledgement

Activity :
 . Select The
 . Enter Com
 . Enter E-Sig

List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0023 Product Name: Aspirin

Specification Number: SNSS 001 AR Number: ASS220068

Batch Number: 44 Packing Type: Box

Storage Conditions: 40°C ± 2°C / 7 H

Qty Withdraw: 2 UOM: Gm

Start Date: 15/07/2022

Used For: Sampling

Comments*: Approved


Validate E-Signature

E-Signature Password* :

LIMS 2.0

6. Enter E-Signature Password

SSNo : PSS-LIMS-SS-006



PHARMA SOFT **STABILITY MANAGEMENT** LIMS APPROVER

Sample Acknowledgement

Activity :
 . Select The
 . Enter Com
 . Enter E-Sig

List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0023 Product Name: Aspirin

Specification Number: SNSS 001 AR Number: ASS220068

Batch Number: 44 Packing Type: Box

Storage Conditions: 40°C ± 2°C / 7 H

Qty Withdraw: 2 UOM: Gm

Start Date: 15/07/2022


Used For: Sampling

Comments*: Approved

Validate E-Signature


E-Signature Password* :

LIMS 2.0

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Sample Acknowledgement
	Test Case Number	LIMS-SM-014
	Actual Result	Pass

7. Click on Submit

SSNo : PSS-LIMS-SS-007



PHARMA SOFT **STABILITY MANAGEMENT** **LIMS APPROVER**

Sample Acknowledgement

Activity :
 . Select The
 . Enter Com
 . Enter E-Sig

List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0023 Product Name: Aspirin

Specification Number: SNSS 001 AR Number: ASS220068

Batch Number: 44 Packing Type: Box

Storage Conditions: 40°C ± 2°C / 75%RH ± 5% UOM: Gm

Qty Withdraw: 2

Start Date: 15/07/2022

Used For: Sampling

Comments*: Approved

Validate E-Signature

E-Signature Password* : *****

Saving...


Submit Cancel

Approve Revert

LIMS 2.0

8. Click on OK button

SSNo : PSS-LIMS-SS-008



PHARMA SOFT **STABILITY MANAGEMENT** **LIMS APPROVER**

Sample Acknowledgement

Activity :
 . Select The
 . Enter Com
 . Enter E-Sig

List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0023 Product Name: Aspirin

Specification Number: SNSS 001 AR Number: ASS220068

Batch Number: 44 Packing Type: Box

Storage Conditions: 40°C ± 2°C / 75%RH ± 5% UOM: Gm

Qty Withdraw: 2

Start Date: 15/07/2022

Used For: Sampling

Comments*: Approved


Success

Submitted successfully.

OK

Approve Revert


LIMS 2.0

	Test Case Name	Approve Sample Acknowledgement
	Test Case Number	LIMS-SM-014
	Actual Result	Pass

9. Click on User Name

SSNo : PSS-LIMS-SS-009

PHARMA SOFT SOL
STABILITY MANAGEMENT


LIMS APPROVER

Sample Acknowledgment Approval

Activity :
 . Select The Record and Click On view.
 . Enter Comments and Click on Approve / Revert.
 . Enter E-Signature Password and Click on Submit.

S.No	Protocol Number	Product Name	Specification Number	AR Number	Storage Condition	Stability Condition	Packing type	Batch Number	Start Date	Done By/On	Logout	Action
1	PS0014	Aspirin	SNSS 001	ASS220038	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	Product Batch No	21/06/2022	LIMS USER / 21/06/2022 14:35:27		View
2	PS0014	Aspirin	SNSS 001	ASS220041	30°C ± 2°C / 65%RH ± 5%RH	INTERMEDIATE	Box	Product Batch No	21/06/2022	LIMS USER / 21/06/2022 14:40:00		View

Go to page:
Row count:
Showing 1-2 of 2

LIMS 2.0
All © rights Reserved PharmaSoftSol

10. Click on Logout

SSNo : PSS-LIMS-SS-010

LIMS



User Name

Password

Module

[Login](#)
[Reset](#)