	Test Case Name	Approve Protocol
	Test Case Number	LIMS-SM-005
	Actual Result	Pass

1.Enter the username, password, Select Module and click on login button


SSNo : PSS-LIMS-SS-001



The login form is titled "LIMS" and features the PharmaSoftSol logo. It contains three input fields: "User Name" with the value "LIMS_APPROVER", "Password" with masked characters "*****", and "Module" with a dropdown menu set to "STABILITY MANAGEMENT". Below these fields are two buttons: "Login" and "Reset".


2.Click on Protocol

SSNo : PSS-LIMS-SS-002



The "Protocol Approval" page shows a table of protocols for approval. The table has columns for S.No., Protocol Number, Product Name, Specification Name, Generic Name, Registered By/On, Request From, and Action. A single record is displayed for protocol PC0014, product PARA65, specification SPPARACHEMICAL, and generic name Paracetamol 650 mg. The "Request From" field shows "LIMS User 31/08/2022 18:24:33" and the "Request From" dropdown is set to "REVERT-RESUBMIT". The "Action" column has a "View & Approval" button. The page also includes a footer with "LIMS 2.0" and "All © rights Reserved PharmaSoftSol".

S.No.	Protocol Number	Product Name	Specification Name	Generic Name	Registered By/On	Request From	Action
1	PC0014	PARA65	SPPARACHEMICAL	Paracetamol 650 mg	LIMS User 31/08/2022 18:24:33	REVERT-RESUBMIT	View & Approval

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Protocol
	Test Case Number	LIMS-SM-005
	Actual Result	Pass

3. Select a Record And Click on View & Approval

SSNo : PSS-LIMS-SS-003

PHARMA SOFT
STABILITY MANAGEMENT
☰
 LIMS Approver

Protocol Objective

Procedure :

↶ ↷
Formats ▾
B
I

Protocol Procedure

Location:

Plant-1

Specification Number:

SPPARACH001

Generic Name:

Paracetamol 650 mg

Test Planner Version:

0

Product Name:

PARA65

Specification Name:

SPPARACHEMICAL

STP Number:

STPCHPARA001

Customer Name:

BA

4. Enter Comments

SSNo : PSS-LIMS-SS-004

PHARMA SOFT
STABILITY MANAGEMENT
☰
 LIMS Approver


Packing Type	Product Batch No	Batch Size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Interval Type
No data available!								

Comments*:

Approved


Approve

Revert

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Protocol
	Test Case Number	LIMS-SM-005
	Actual Result	Pass

5. Click on Approve


SSNo : PSS-LIMS-SS-005




The screenshot shows the PHARMA SOFT STABILITY MANAGEMENT interface. At the top, there is a header with the PHARMA SOFT logo and the text 'STABILITY MANAGEMENT'. Below the header, there is a table with columns: Packing Type, Product Batch No, Batch Size, Manufactured Date, Expiry Date, Storage Condition, Orientation, Stability Condition, and Interval Type. The table is currently empty, with a message 'No data available!' displayed below it. In the center of the screen, there is a 'Validate E-Signature' dialog box. The dialog box has a title bar 'Validate E-Signature' and a label 'E-Signature Password*:' followed by a text input field. Below the input field, there are two buttons: 'Submit' and 'Cancel'. At the bottom of the main window, there are two buttons: 'Approve' (green) and 'Revert' (red). The 'Comments*' field contains the text 'Approved'.

6. Enter E-Signature Password

SSNo : PSS-LIMS-SS-006



This screenshot is similar to the previous one, showing the same PHARMA SOFT STABILITY MANAGEMENT interface. The 'Validate E-Signature' dialog box is still open, but the 'E-Signature Password*' input field now contains several asterisks, indicating that a password has been entered. The 'Approve' and 'Revert' buttons remain at the bottom of the main window, and the 'Comments*' field still shows 'Approved'.

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Approve Protocol
	Test Case Number	LIMS-SM-005
	Actual Result	Pass

7.Click on Submit

SSNo : PSS-LIMS-SS-007

PHARMA SOFT
STABILITY MANAGEMENT
☰
👤 LIMS Approver ▾

Packing Type	Product Batch No	Batch Size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Interval Type
No data available!								

Comments*:

Approved

Approve
Revert

8.Click on OK button

SSNo : PSS-LIMS-SS-008

PHARMA SOFT
STABILITY MANAGEMENT
☰
👤 LIMS Approver ▾

Storage Condition Interval Details


Packing Type	Product Batch No	Batch Size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Interval Type	Int
No data available!									

Comments*:

Approved

Success
Submitted successfully.
OK


Approve
Revert

	Test Case Name	Approve Protocol
	Test Case Number	LIMS-SM-005
	Actual Result	Pass

9.Click on User Name


SSNo : PSS-LIMS-SS-009

PHARMA SOFT SOL
STABILITY MANAGEMENT

 LIMS Approver

Protocol Approval

Activity :
 . Select The Record and Click On View & Approval.
 . Enter Comments and Click on Approve / Revert.
 . Enter E-Signature Password and Click on Submit.



List of Protocols for Approval

S.No.	Protocol Number	Product Name	Specification Name	Generic Name	Registered By/On	Request From	Action
No data available!							

Row count: 10

LIMS 2.0
 All © rights Reserved PharmaSoftSol

10.Click on Logout

SSNo : PSS-LIMS-SS-010

LIMS



User Name

Password

Module

Select

Login

Reset