
	<b>Test Case Name</b>	Create Schedule Test
	<b>Test Case Number</b>	LIMS-SM-007
	<b>Actual Result</b>	Pass

1. Enter the username, password, Select Module and click on login button


SSNo : PSS-LIMS-SS-001




The image shows the LIMS login interface. It features the PharmaSoftSol logo at the top left. The main heading is 'LIMS'. Below it, there are three input fields: 'User Name' with the value 'LIMS\_USER', 'Password' with masked characters '\*\*\*\*\*', and 'Module' with a dropdown menu set to 'STABILITY MANAGEMENT'. At the bottom of the form are two buttons: 'Login' and 'Reset'.

2. Click on Schedule Test

SSNo : PSS-LIMS-SS-002



The screenshot displays the 'Schedule Test Selection' page within the 'STABILITY MANAGEMENT' module. The page has a header with 'PHARMA SOFT' and 'STABILITY MANAGEMENT' logos, and a user profile 'LIMS User'. The main content area includes a 'Schedule Test Selection' tab, an 'Activity' section with instructions, and a 'List of Reverted Schedules' table. The table has columns for S.No., Protocol Number, Product Name, Generic Name, Specification Name, Reverted By, Reverted On, and Comments. Below the table, it indicates 'No data available!'. The footer shows 'LIMS 2.0' and 'All © rights Reserved PharmaSoftSol'.

	<b>Test Case Name</b>	Create Schedule Test
	<b>Test Case Number</b>	LIMS-SM-007
	<b>Actual Result</b>	Pass

3.Click on Create

SSNo : PSS-LIMS-SS-003

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS User

### Schedule Test Selection

Action
Schedule Intervals
Schedule Test Selection

**Activity :**

- . Select Create, Select The Record and Click on Next.
- . If Rejected is Required, Then Select Rejected.

☒ Schedule Test Selection
 ☐ Reverted Schedule Test Selection

#### Schedule Create

S.No.	Protocol Number	Product	Specification Name	Generic Name
No data available!				

Go to page: 1
Row count: 10

IMS 2.0
 All © rights Reserved PharmaSoftSol

4.Select a Record

SSNo : PSS-LIMS-SS-004

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS User

### Schedule Test Selection

Action
Schedule Intervals
Schedule Test Selection

**Activity :**

- . Select Create, Select The Record and Click on Next.
- . If Rejected is Required, Then Select Rejected.

☒ Schedule Test Selection
 ☐ Reverted Schedule Test Selection


#### Schedule Create

S.No.	Protocol Number	Product	Specification Name	Generic Name
1	PC0011	PARA65	SPPARACHEMICAL	Paracetamol 650 mg

<<
 <
 1
 >
 >>

Go to page: 1
Row count: 10

Showing 1-1 of 1

	<b>Test Case Name</b>	Create Schedule Test
	<b>Test Case Number</b>	LIMS-SM-007
	<b>Actual Result</b>	Pass

5. Click on Next

SSNo : PSS-LIMS-SS-005

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS User

### Schedule Test Selection

Action
Schedule Intervals
Schedule Test Selection

Record and Click on Next.  
When Select Rejected.

☐
Reverted Schedule Test Selection

Product	Specification Name	Generic Name
PARA65	SPPARACHEMICAL	Paracetamol 650 mg

Page: 1 Row count: 10 Showing 1-1 of 1

. Select The Record and Click On Next.

**Schedule Test Selection**

☐ S.No ☐ Packing Type

Go to page: 1 Row count: 10

LIMS 2.0 All © rights Reserved PharmaSoftSol

6. Select a Record

SSNo : PSS-LIMS-SS-006

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS User

### Schedule Test Selection


Action
Schedule Intervals
Schedule Test Selection

. Select The Record and Click On Next.

S.No	Packing Type	Product Batch No	Mfg.Date	Exp.Date	Storage Conditions	Orientation	Stability Conditions	Schedule Interval Type
1	Blister	44	02-11-2020	02-11-2021	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH

<< < 1 > >> Go to page: 1 Row count: 10 Showing 1-1 of 1

LIMS 2.0 All © rights Reserved PharmaSoftSol

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Create Schedule Test
	<b>Test Case Number</b>	LIMS-SM-007
	<b>Actual Result</b>	Pass

7. Click on Next

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS User

### Schedule Test Selection

Action
Schedule Intervals
Schedule Test Selection

. Select The Test, Select Approval Life Cycle and Click On Finish.  
. Enter E-Signature Password and Click On Submit.

Set1

Packing Type: Blister
Storage Conditions: 40°C ± 2°C / 75%RH ± 5%RH
Stability Condition: ACCELERATED
Product Batch No: 44
Schedule Interval Type: MONTH
Orientation: UP RIGHT

**Interval and Test Details**

S.No	Test Name	
1	PARA65CHDescription	<input checked="" type="checkbox"/>

IMS 2.0
All © rights Reserved PharmaSoftSol

8. Select Approval Cycle

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS User

### Schedule Test Selection

Action
Schedule Intervals
Schedule Test Selection

. Select The Test, Select Approval Life Cycle and Click On Finish.  
. Enter E-Signature Password and Click On Submit.

Set1


Packing Type: Blister
Storage Conditions: 40°C ± 2°C / 75%RH ± 5%RH
Stability Condition: ACCELERATED
Product Batch No: 44
Schedule Interval Type: MONTH
Orientation: UP RIGHT

**Interval and Test Details**

S.No	Test Name	
1	PARA65CHDescription	<input checked="" type="checkbox"/>
1	Number of Syringes required at Each Point	50


Approval Life Cycle\*: ApproverLC View

IMS 2.0
All © rights Reserved PharmaSoftSol

	<b>Test Case Name</b>	Create Schedule Test
	<b>Test Case Number</b>	LIMS-SM-007
	<b>Actual Result</b>	Pass


9. Click on Finish

SSNo : PSS-LIMS-SS-009




10. Enter E-Signature Password

SSNo : PSS-LIMS-SS-010



SSNo : PSS-LIMS-SS-010

Page 6 of 8

	<b>Test Case Name</b>	Create Schedule Test
	<b>Test Case Number</b>	LIMS-SM-007
	<b>Actual Result</b>	Pass

11.Click on OK button

SSNo : PSS-LIMS-SS-011

**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS User

Packing Type: Blister Product Batch No: 44  
Storage Conditions: 40°C ± 2°C / 75%RH ± 5%RH Schedule Interval Type: MONTH  
Stability Condition: ACCELERATED Orientation: UP RIGHT

**Interval and Test Details**

S.No	Test Name
1	PARA65CHDescription

Success  
Submitted successfully.  
OK

Approval Life Cycle\*: ApproverLC

Previous Finish

LIMS 2.0 All © rights Reserved PharmaSoftSol

12.Click on User Name

SSNo : PSS-LIMS-SS-012

**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS User

☐ Schedule Test Selection ☒ Reverted Schedule Test Selection

**List of Reverted Schedules**


S.No.	Protocol Number	Product Name	Generic Name	Specification Name	Reverted By	Reverted On	Comments
No data available!							

Row count: 10

Logout


Previous Next

LIMS 2.0 All © rights Reserved PharmaSoftSol


	<b>Test Case Name</b>	Create Schedule Test
	<b>Test Case Number</b>	LIMS-SM-007
	<b>Actual Result</b>	Pass

13. Click on Logout

SSNo : PSS-LIMS-SS-013



**LIMS**



User Name

Password

Module