
 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

1. Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001




The screenshot shows the LIMS login interface. At the top center, the word "LIMS" is displayed in yellow. Below it is the PharmaSoftSol logo. The login form contains the following fields:

- User Name:** A text input field containing "LIMS_USER".
- Password:** A password input field with masked characters "*****".
- Module:** A dropdown menu currently showing "STABILITY MANAGEMENT".
- Login/Reset Buttons:** Two blue buttons labeled "Login" and "Reset" are positioned at the bottom of the form.

2. Click on Loading

SSNo : PSS-LIMS-SS-002



The screenshot displays the "Loading" section of the LIMS application. The header includes the PharmaSoftSol logo, the text "STABILITY MANAGEMENT", and a user profile icon labeled "LIMS USER".

Loading Section:

- Activity:**
 - . Select Create as Loading, Select The Record and Click on Next.
 - . If Reverted Loading by Approver is Required, Then Select Reverted Loading by Approver.
- Radio Buttons:**
 - ☒ Loading
 - ☐ Reverted Loading by Approver


Loading Details Table:

Protocol Number	Product	Generic Name	Specification Name	Protocol
PS0021	Aspirin	Aspirin and codeine	SS Aspirin Spec	View
PS0022	Aspirin	Aspirin and codeine	SS Aspirin Spec	View
PS0023	Aspirin	Aspirin and codeine	SS Aspirin Spec	View

At the bottom of the table, there is a pagination bar with the following elements:

- Navigation buttons: <<, <, 1, >, >>
- Go to page: 1
- Row count: 10
- Showing 1-3 of 3

The footer of the application shows "IMS 2.0" on the left and "All © rights Reserved PharmaSoftSol" on the right.

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

3. Click on Rejected

SSNo : PSS-LIMS-SS-003

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading

Activity :
. Select The Record and Click on Next.

☐ Loading
 ☒ Reverted Loading by Approver

Protocol Number	Product Name	Generic Name	Specification Name	Reverted By/On	Comments
PS0024	Aspirin	Aspirin and codeine	SS Aspirin Spec	LIMS APPROVER 15/07/2022 11:08:43	View

<< < 1 > >>
Go to page: 1
Row count: 10
Showing 1-1 of 1

IMS 2.0
All © rights Reserved PharmaSoftSol

4. Select a Record

SSNo : PSS-LIMS-SS-004

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading


Activity :
. Select The Record and Click on Next.

☐ Loading
 ☒ Reverted Loading by Approver

Protocol Number	Product Name	Generic Name	Specification Name	Reverted By/On	Comments
PS0024	Aspirin	Aspirin and codeine	SS Aspirin Spec	LIMS APPROVER 15/07/2022 11:08:43	View

<< < 1 > >>
Go to page: 1
Row count: 10
Showing 1-1 of 1

IMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

5. Click on Next

SSNo : PSS-LIMS-SS-005

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action Loading

Reverted By/On

Comments

LIMS APPROVER 15/07/2022 11:08:43

View

Showing 1-1 of 1

. Modify The Required Fields and Click On Finish.

. Enter E-Signature Password and Click On Submit.

Loading Details

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Rel Coi
No data available!									

IMS 2.0

All © rights Reserved PharmaSoftSol

6. Clear Analysis Quantity

SSNo : PSS-LIMS-SS-006

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action Loading

. Modify The Required Fields and Click On Finish.

. Enter E-Signature Password and Click On Submit.


Loading Details

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH		3	Gm	3	2	Gm	Evaluate	15-07-2022	Select

Approval From*: LIMS APPROVER Select

IMS 2.0

All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

7. Enter Analysis Quantity

SSNo : PSS-LIMS-SS-007

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading

. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	3	Gm	3	2	Gm	Evaluate	15-07-2022	Select

Approval From*:
LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

8. Clear Reference Quantity

SSNo : PSS-LIMS-SS-008

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading


Action
Loading

. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4		Gm	3	2	Gm	Evaluate	15-07-2022	Select

Approval From*:
LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

9. Enter Reference Quantity

SSNo : PSS-LIMS-SS-009

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading

. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	3	2	Gm	Evaluate	15-07-2022	Select

Approval From*:
LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

10. Select UOM

SSNo : PSS-LIMS-SS-010

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading


Action
Loading

. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	3	2	Gm	Evaluate	15-07-2022	Select

Approval From*:
LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

11. Clear Analysis Quantity

SSNo : PSS-LIMS-SS-011

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading

. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm		2	Gm	Evaluate	15-07-2022	Select

Approval From*:

LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

12. Enter Analysis Quantity Count

SSNo : PSS-LIMS-SS-012

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading


. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2	2	Gm	Evaluate	15-07-2022	Select

Approval From*:

LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

13. Clear Reference Quantity

SSNo : PSS-LIMS-SS-013

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading

. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2		Gm	Evaluate	15-07-2022	Select

Approval From*:

LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

14. Enter Reference Quantity Count

SSNo : PSS-LIMS-SS-014

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading


. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2		Gm	Evaluate	15-07-2022	Select

Approval From*:

LIMS APPROVER
Select

IMS 2.0
All © rights Reserved PharmaSoftSol

	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

15.Select UOM

SSNo : PSS-LIMS-SS-015

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading

Action
Loading

. Modify The Required Fields and Click On Finish.
. Enter E-Signature Password and Click On Submit.

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	Evaluate	15-07-2022	Select

Approval From*: LIMS APPROVER Select

IMS 2.0
All © rights Reserved PharmaSoftSol

16.Click on Finish

SSNo : PSS-LIMS-SS-016

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading Details

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	Evaluate	15-07-2022	Select


Approval From*: LIMS APPROVER

Validate E-Signature
E-Signature Password*:
Submit
Cancel

Previous

Finish

IMS 2.0
All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

17. Enter E-Signature Password

SSNo : PSS-LIMS-SS-017

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	Evaluate	15-07-2022	Select

Approval From*: LIMS APPROVER

Validate E-Signature

E-Signature Password*:

Submit Cancel

Previous Finish

LIMS 2.0 All © rights Reserved PharmaSoftSol

18. Click on Submit

SSNo : PSS-LIMS-SS-018

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	Evaluate	15-07-2022	Select

Approval From*: LIMS APPROVER

Validate E-Signature


E-Signature Password*:

Saving...

Submit Cancel

Previous Finish

LIMS 2.0 All © rights Reserved PharmaSoftSol

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

19. Click on OK button

SSNo : PSS-LIMS-SS-019

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Loading Details

Packing Type	Batch Number	Storage Conditions	Stability Conditions	Schedule Interval Type	Analysis Quantity	Reference Quantity	UOM (Quantity)	Analysis Qty Count	Ref. Quantity Count	UOM (Qty Count)	Total Count	Start Date	Chamber/Rack
Box	44	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	MONTH	4	1	Gm	2	3	Gm	Evaluate	15-07-2022	Select

Approval From*: LIMS APPROVER

Success

Submitted successfully.

OK

Previous Finish

LIMS 2.0
All © rights Reserved PharmaSoftSol

20. Click on User Name

SSNo : PSS-LIMS-SS-020

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

☒ Loading
☐ Reverted Loading by Approver

Loading Details


Protocol Number	Product	Generic Name	Specification Name	Protocol
PS0021	Aspirin	Aspirin and codeine	SS Aspirin Spec	View
PS0022	Aspirin	Aspirin and codeine	SS Aspirin Spec	View
PS0023	Aspirin	Aspirin and codeine	SS Aspirin Spec	View

[<<](#)
[<](#)
[1](#)
[>](#)
[>>](#)
Go to page Row count:

Showing 1-3 of 3 Logout


Previous Next

LIMS 2.0
All © rights Reserved PharmaSoftSol


	Test Case Name	Reinitiate Rejected Loading
	Test Case Number	PSS-LIMS-019
	Actual Result	Pass

21.Click on Logout

SSNo : PSS-LIMS-SS-021



LIMS



User Name

Password

Module