	<b>Test Case Name</b>	Reject Loading
	<b>Test Case Number</b>	PSS-LIMS-018
	<b>Actual Result</b>	Pass

1. Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001



The screenshot shows the LIMS login interface. At the top center, the word "LIMS" is displayed in yellow. Below it is the PharmaSoftSol logo. The login form contains three input fields: "User Name" with the value "LIMS\_APPROVER", "Password" with masked characters "\*\*\*\*\*", and "Module" with a dropdown menu set to "STABILITY MANAGEMENT". At the bottom of the form are two buttons: "Login" and "Reset".

1. Click on Loading

SSNo : PSS-LIMS-SS-001



The screenshot shows the "Loading Approval" screen in the LIMS system. The header includes the PharmaSoftSol logo, the text "STABILITY MANAGEMENT", and a user profile icon labeled "LIMS Approver". Below the header, there is a section titled "Loading Approval" with instructions: "Activity : . Select The Record and Click On View. . Enter Comments and Click on Approve / Revert. . Enter E-Signature Password and Click on Submit." Below this is a table titled "Loading Details" with columns: S.No., Protocol Number, Product Name, Generic Name, Specification Name, Loading By/On, and Details. The table contains one row with the following data: S.No. 1, Protocol Number PC0015, Product Name PARA65, Generic Name Paracetamol 650 mg, Specification Name SPPARACHEMICAL, Loading By/On LIMS User 31/08/2022 20:42:53, and a "View" button. Below the table is a pagination bar with navigation arrows, "Go to page 1", "Row count: 10", and "Showing 1-1 of 1". At the bottom of the screen, there is a footer with "LIMS 2.0" and "All © rights Reserved PharmaSoftSol".

2. Select a Record And Click on View

SSNo : PSS-LIMS-SS-002

<b>Test Case Name</b>	Reject Loading
<b>Test Case Number</b>	PSS-LIMS-018
<b>Actual Result</b>	Pass

### Loading Details

**Activity :**


- Select The Record as per the Requirement.
- Enter Comments and Remarks.
- Enter E-Signature Panel.

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	2	3	Gm	3	2	Gm	31/08/2022	<button style="background-color: #0056b3; color: white; padding: 2px 5px;">View</button>

**Comments\*:**

Maximum allowed characters should be 1500

Approve
Revert

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Loading
	<b>Test Case Number</b>	PSS-LIMS-018
	<b>Actual Result</b>	Pass

### 3. Enter Comments

SSNo : PSS-LIMS-SS-003

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver

Loading Approval

Activity :  
. Select The Record at  
. Enter Comments and  
. Enter E-Signature Pa

**Loading Details**  
S.No. Protocol  
1 PC0015  
<< < 1 > >> Go to p

**Loading Details**

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	2	3	Gm	3	2	Gm	31/08/2022	<a href="#">View</a>

Comments\* :

Rejected

Approve
Revert

Details  
View  
Showing 1-1 of 1

LIMS 2.0
All © rights Reserved PharmaSoftSol

### 4. Click on Reject

SSNo : PSS-LIMS-SS-004

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver

Loading Approval

Activity :  
. Select The Record at  
. Enter Comments and  
. Enter E-Signature Pa

**Loading Details**  
S.No. Protocol  
1 PC0015  
<< < 1 > >> Go to p

**Loading Details**

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Analysis Qty Count	Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP R								Gm	31/08/2022	<a href="#">View</a>

Comments\* :

Rejected


Validate E-Signature

E-Signature Password\* :

Submit
Cancel

Details  
View  
Showing 1-1 of 1

LIMS 2.0
All © rights Reserved PharmaSoftSol

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Loading
	<b>Test Case Number</b>	PSS-LIMS-018
	<b>Actual Result</b>	Pass

## 5. Enter E-Signature Password

SSNo : PSS-LIMS-SS-005

**PHARMA SOFT** **STABILITY MANAGEMENT** ☰ 👤 LIMS Approver

**Loading Approval**

Activity :  
 . Select The Record at  
 . Enter Comments and  
 . Enter E-Signature Pa

**Loading Details**

S.No. Protocol  
 1 PC0015

<< < 1 > >> Go to p

**Loading Details**

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP R								Gm	31/08/2022	<a href="#">View</a>

Comments \* : Rejected

**Validate E-Signature**

E-Signature Password\* :

[Submit](#) [Cancel](#)

**LIMS 2.0** All © rights Reserved PharmaSoftSol

## 6. Click on Submit

SSNo : PSS-LIMS-SS-006

**PHARMA SOFT** **STABILITY MANAGEMENT** ☰ 👤 LIMS Approver

**Loading Approval**

Activity :  
 . Select The Record at  
 . Enter Comments and  
 . Enter E-Signature Pa

**Loading Details**

S.No. Protocol  
 1 PC0015

<< < 1 > >> Go to p

**Loading Details**

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/f Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP R								Gm	31/08/2022	<a href="#">View</a>


Comments \* : Rejected

**Validate E-Signature**

E-Signature Password\* :

[Saving...](#) [Submit](#) [Cancel](#)

**LIMS 2.0** All © rights Reserved PharmaSoftSol

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Loading
	<b>Test Case Number</b>	PSS-LIMS-018
	<b>Actual Result</b>	Pass

7. Click on OK button

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver

Loading Approval

Activity :  
. Select The Record and Click On View.  
. Enter Comments and Click on Approve / Revert.  
. Enter E-Signature Password and Click on Submit.

Loading Details
S.No. Protocol  
1 PC0015

Loading Details

Packing Type	Product Batch No	Storage Conditions	orientation	Stability Conditions	Schedule Interval	Analysis Quantity	Reference Quantity	UOM(Quantity)	Qty Count	Analysis Reference Qty Count	UOM(Qty Count)	Start Date	Chamber/F Details
Blister	44	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	2	3	Gm	3	2	Gm	31/08/2022	<a href="#">View</a>

Comments\* : Rejected

Success

Submitted successfully.

OK

Approve

Revert

LIMS 2.0
All © rights Reserved PharmaSoftSol

8. Click on User Name

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver


Loading Approval

Activity :  
. Select The Record and Click On View.  
. Enter Comments and Click on Approve / Revert.  
. Enter E-Signature Password and Click on Submit.

Loading Details
S.No. Protocol Number Product Name Generic Name Specification Name Loading By/On Logout Details


Row count: 10

LIMS 2.0
All © rights Reserved PharmaSoftSol

	<b>Test Case Name</b>	Reject Loading
	<b>Test Case Number</b>	PSS-LIMS-018
	<b>Actual Result</b>	Pass

9.Click on Logout  
SSNo : PSS-LIMS-SS-009

LIMS



User Name

User Name

Password

Password

Module

Select

Login

Reset