


|   |                         |                 |
|---|-------------------------|-----------------|
|  | <b>Test Case Name</b>   | Approve Loading |
|   | <b>Test Case Number</b> | LIMS-SM-012     |
|   | <b>Actual Result</b>    | Pass            |

1.Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001



1.Click on Loading


SSNo : PSS-LIMS-SS-001



2.Select a Record And Click on View

SSNo : PSS-LIMS-SS-002




|  |                         |                 |
|--|-------------------------|-----------------|
|  <b>PHARMA<br/>SOFT SOL</b><br><small>COMMITTED TO PAPERLESS COMPLIANCE</small> | <b>Test Case Name</b>   | Approve Loading |
|  | <b>Test Case Number</b> | LIMS-SM-012     |
|  | <b>Actual Result</b>    | Pass            |

### 3. Enter Comments

SSNo : PSS-LIMS-SS-003

**PHARMA SOFT**

**STABILITY MANAGEMENT**

 LIMS Approver

Loading Approval

Activity :  
. Select The Record at  
. Enter Comments and  
. Enter E-Signature Pa

**Loading Details**  
S.No. Protocol  
1 PC0011  
<< < 1 > >> Go to p

Loading Details

| Packing Type | Product Batch No | Storage Conditions        | orientation | Stability Conditions | Schedule Interval | Analysis Quantity | Reference Quantity | UOM(Quantity) | Analysis Qty Count | Reference Qty Count | UOM(Qty Count) | Start Date | Chamber/f Details    |
|--------------|------------------|---------------------------|-------------|----------------------|-------------------|-------------------|--------------------|---------------|--------------------|---------------------|----------------|------------|----------------------|
| Blister      | 44               | 40°C ± 2°C / 75%RH ± 5%RH | UP RIGHT    | ACCELERATED          | MONTH             | 2                 | 3                  | Gm            | 3                  | 2                   | Gm             | 31/08/2022 | <a href="#">View</a> |

Comments \*

Approved

[Approve](#)
[Revert](#)

Details  
6:18:01  
[View](#)  
Showing 1-1 of 1


LIMS 2.0
All © rights Reserved PharmaSoftSol

### 4. Click on Approve

SSNo : PSS-LIMS-SS-004

**PHARMA SOFT**

**STABILITY MANAGEMENT**

 LIMS Approver

Loading Approval

Activity :  
. Select The Record at  
. Enter Comments and  
. Enter E-Signature Pa

**Loading Details**  
S.No. Protocol  
1 PC0011  
<< < 1 > >> Go to p

Loading Details

| Packing Type | Product Batch No | Storage Conditions        | orientation | Stability Conditions | Schedule Interval | Analysis Quantity | Reference Quantity | UOM(Quantity) | Analysis Qty Count | Reference Qty Count | UOM(Qty Count) | Start Date | Chamber/f Details    |
|--------------|------------------|---------------------------|-------------|----------------------|-------------------|-------------------|--------------------|---------------|--------------------|---------------------|----------------|------------|----------------------|
| Blister      | 44               | 40°C ± 2°C / 75%RH ± 5%RH | UP R        |                      |                   |                   |                    |               |                    |                     | Gm             | 31/08/2022 | <a href="#">View</a> |

Comments \*

Appro


Validate E-Signature

E-Signature Password\* :

Submit Cancel

Details  
6:18:01  
[View](#)  
Showing 1-1 of 1

LIMS 2.0
All © rights Reserved PharmaSoftSol

|  |                         |                 |
|--|-------------------------|-----------------|
|  <b>PHARMA<br/>SOFT SOL</b><br><small>COMMITTED TO PAPERLESS COMPLIANCE</small> | <b>Test Case Name</b>   | Approve Loading |
|  | <b>Test Case Number</b> | LIMS-SM-012     |
|  | <b>Actual Result</b>    | Pass            |

## 5. Enter E-Signature Password

SSNo : PSS-LIMS-SS-005

**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS Approver

**Loading Approval**

Activity :  
 . Select The Record at  
 . Enter Comments and  
 . Enter E-Signature Pa

**Loading Details**

S.No. Protocol  
 1 PC0011

<< < 1 > >> Go to p

**Loading Details**

| Packing Type | Product Batch No | Storage Conditions        | orientation | Stability Conditions | Schedule Interval | Analysis Quantity | Reference Quantity | UOM(Quantity) | Qty Count | Reference Qty Count | UOM(Qty Count) | Start Date | Chamber/f Details    |
|--------------|------------------|---------------------------|-------------|----------------------|-------------------|-------------------|--------------------|---------------|-----------|---------------------|----------------|------------|----------------------|
| Blister      | 44               | 40°C ± 2°C / 75%RH ± 5%RH | UP R        |                      |                   |                   |                    |               |           |                     | Gm             | 31/08/2022 | <a href="#">View</a> |

Comments\* : Appro

**Validate E-Signature**

E-Signature Password\* :

[Submit](#) [Cancel](#)

LIMS 2.0 All © rights Reserved PharmaSoftSol

## 6. Click on Submit

SSNo : PSS-LIMS-SS-006

**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS Approver

**Loading Approval**

Activity :  
 . Select The Record at  
 . Enter Comments and  
 . Enter E-Signature Pa

**Loading Details**

S.No. Protocol  
 1 PC0011

<< < 1 > >> Go to p

**Loading Details**

| Packing Type | Product Batch No | Storage Conditions        | orientation | Stability Conditions | Schedule Interval | Analysis Quantity | Reference Quantity | UOM(Quantity) | Qty Count | Reference Qty Count | UOM(Qty Count) | Start Date | Chamber/f Details    |
|--------------|------------------|---------------------------|-------------|----------------------|-------------------|-------------------|--------------------|---------------|-----------|---------------------|----------------|------------|----------------------|
| Blister      | 44               | 40°C ± 2°C / 75%RH ± 5%RH | UP R        |                      |                   |                   |                    |               |           |                     | Gm             | 31/08/2022 | <a href="#">View</a> |


Comments\* : Appro

**Validate E-Signature**

E-Signature Password\* :

[Saving...](#)  
[Submit](#) [Cancel](#)

LIMS 2.0 All © rights Reserved PharmaSoftSol

|  |                         |                 |
|--|-------------------------|-----------------|
|  <b>PHARMA<br/>SOFT SOL</b><br><small>COMMITTED TO PAPERLESS COMPLIANCE</small> | <b>Test Case Name</b>   | Approve Loading |
|  | <b>Test Case Number</b> | LIMS-SM-012     |
|  | <b>Actual Result</b>    | Pass            |

7. Click on OK button

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver

Loading Approval

Activity :  
. Select The Record and Click On View.  
. Enter Comments and Click on Approve / Revert.  
. Enter E-Signature Password and Click on Submit.

Loading Details
S.No. Protocol  
1 PC0011

Loading Details

| Packing Type | Product Batch No | Storage Conditions        | orientation | Stability Conditions | Schedule Interval | Analysis Quantity | Reference Quantity | UOM(Quantity) | Qty Count | Analysis Reference Qty Count | UOM(Qty Count) | Start Date | Chamber/F Details    |
|--------------|------------------|---------------------------|-------------|----------------------|-------------------|-------------------|--------------------|---------------|-----------|------------------------------|----------------|------------|----------------------|
| Blister      | 44               | 40°C ± 2°C / 75%RH ± 5%RH | UP RIGHT    | ACCELERATED          | MONTH             | 2                 | 3                  | Gm            | 3         | 2                            | Gm             | 31/08/2022 | <a href="#">View</a> |

Comments : Approved

Success

Submitted successfully.

OK

[Approve](#)
[Revert](#)

6:18:01
[View](#)

Showing 1-1 of 1

LIMS 2.0
All © rights Reserved PharmaSoftSol

8. Click on User Name

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS Approver

Loading Approval


Activity :  
. Select The Record and Click On View.  
. Enter Comments and Click on Approve / Revert.  
. Enter E-Signature Password and Click on Submit.

Loading Details
S.No. Protocol Number Product Name Generic Name Specification Name Loading By/On Logout Details

Row count: 10


No data available!

LIMS 2.0
All © rights Reserved PharmaSoftSol

|   |                         |                 |
|---|-------------------------|-----------------|
|  | <b>Test Case Name</b>   | Approve Loading |
|   | <b>Test Case Number</b> | LIMS-SM-012     |
|   | <b>Actual Result</b>    | Pass            |

9.Click on Logout  
SSNo : PSS-LIMS-SS-009

LIMS



User Name

User Name

Password

Password

Module

Select

LoginReset