	<b>Test Case Name</b>	Approve Schedule Test
	<b>Test Case Number</b>	LIMS-SM-008
	<b>Actual Result</b>	Pass

1.Enter the username, password, Select Module and click on login button


SSNo : PSS-LIMS-SS-001



The screenshot shows the LIMS login interface. At the top, the 'LIMS' title is displayed in yellow. Below it, the PharmaSoftSol logo is visible. The login form contains three fields: 'User Name' with the value 'LIMS\_APPROVER', 'Password' with masked characters '\*\*\*\*\*', and 'Module' with a dropdown menu set to 'STABILITY MANAGEMENT'. At the bottom of the form are two buttons: 'Login' and 'Reset'.


2.Click on Protocol

SSNo : PSS-LIMS-SS-002



The screenshot displays the 'Schedule Approval' screen. At the top, the PharmaSoftSol logo and 'STABILITY MANAGEMENT' title are shown. The user 'LIMS APPROVER' is logged in. The screen includes instructions for the activity: 'Select The Record', 'Enter Comments and Click on Approve / Revert', and 'Enter E-Signature Password and Click on Submit'. Below this is a table titled 'List of Schedules for Approval' with columns: S.No., Protocol Number, Product Name, Generic Name, Specification Name, Customer, Registered By/On, Schedule Tests, Protocol, and Protocol Modify Details. The table contains two rows of data for Aspirin products. At the bottom, there is a 'Comments' section with a text area and a note about character limits, followed by 'Approve' and 'Revert' buttons.

S.No.	Protocol Number	Product Name	Generic Name	Specification Name	Customer	Registered By/On	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	SS Aspirin Spec	Dr Reddys	LIMS USER 27/06/2022 16:51:02	<a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>
2	PS0023	Aspirin	Aspirin and codeine	SS Aspirin Spec	Dr Reddys	LIMS USER 15/07/2022 10:43:18	<a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Approve Schedule Test
	<b>Test Case Number</b>	LIMS-SM-008
	<b>Actual Result</b>	Pass

### 3.Select a Record

SSNo : PSS-LIMS-SS-003

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
 LIMS APPROVER

#### Schedule Approval

**Activity :**  
 . Select The Record.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.

S.No.	Protocol Number	Product Name	Generic Name	Specification Name	Customer	Registered By/On	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	SS Aspirin Spec	Dr Reddys	LIMS USER 27/06/2022 16:51:02	<a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>
2	PS0023	Aspirin	Aspirin and codeine	SS Aspirin Spec	Dr Reddys	LIMS USER 15/07/2022 10:43:18	<a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>

[<<](#)
[<](#)
[1](#)
[>](#)
[>>](#)
 Go to page:  Row count: 
Showing 1-2 of 2

**Comments \*\*:**

[Approve](#)
[Revert](#)

**LIMS 2.0**
 All © rights Reserved PharmaSoftSol

### 4.Enter Comments

SSNo : PSS-LIMS-SS-004

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
 LIMS APPROVER

#### Schedule Approval

**Activity :**  
 . Select The Record.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.


S.No.	Protocol Number	Product Name	Generic Name	Specification Name	Customer	Registered By/On	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	SS Aspirin Spec	Dr Reddys	LIMS USER 27/06/2022 16:51:02	<a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>
2	PS0023	Aspirin	Aspirin and codeine	SS Aspirin Spec	Dr Reddys	LIMS USER 15/07/2022 10:43:18	<a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>

[<<](#)
[<](#)
[1](#)
[>](#)
[>>](#)
 Go to page:  Row count: 
Showing 1-2 of 2

**Comments \*\*:**

[Approve](#)
[Revert](#)

**LIMS 2.0**
 All © rights Reserved PharmaSoftSol

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Approve Schedule Test
	<b>Test Case Number</b>	LIMS-SM-008
	<b>Actual Result</b>	Pass

5. Click on Approve

SSNo : PSS-LIMS-SS-005

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
 LIMS APPROVER

### Schedule Approval

**Activity :**  
 . Select The Record.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.

S.No.	Protocol Number	Product Name	Generic Name	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	1:02 <a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>
2	PS0023	Aspirin	Aspirin and codeine	8:18 <a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>

<< < 1 > >> Go to page: 1 Row count: 10 Showing 1-2 of 2

**Comments \*:** Approved

[Approve](#) [Revert](#)

**Validate E-Signature**

E-Signature Password\* :

[Submit](#) [Cancel](#)

**LIMS 2.0**
All © rights Reserved PharmaSoftSol

6. Enter E-Signature Password

SSNo : PSS-LIMS-SS-006

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
 LIMS APPROVER

### Schedule Approval

**Activity :**  
 . Select The Record.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.

S.No.	Protocol Number	Product Name	Generic Name	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	1:02 <a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>
2	PS0023	Aspirin	Aspirin and codeine	8:18 <a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>

<< < 1 > >> Go to page: 1 Row count: 10 Showing 1-2 of 2

**Comments \*:** Approved


[Approve](#) [Revert](#)

**Validate E-Signature**

E-Signature Password\* :

[Submit](#) [Cancel](#)

**LIMS 2.0**
All © rights Reserved PharmaSoftSol

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Approve Schedule Test
	<b>Test Case Number</b>	LIMS-SM-008
	<b>Actual Result</b>	Pass

7. Click on Submit

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
 LIMS APPROVER

### Schedule Approval

**Activity :**  
 . Select The Record.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.

S.No.	Protocol Number	Product Name	Generic Name	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	16:02	<a href="#">View</a>	<a href="#">View</a>
2	PS0023	Aspirin	Aspirin and codeine	18:18	<a href="#">View</a>	<a href="#">View</a>

Go to page: 1 Row count: 10

**Comments \***: Approved

[Approve](#)
[Revert](#)

**Validate E-Signature**

E-Signature Password\* :

[Saving...](#)

[Submit](#)
[Cancel](#)

LIMS 2.0
 All © rights Reserved PharmaSoftSol

8. Click on OK button

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
 LIMS APPROVER

### Schedule Approval

**Activity :**  
 . Select The Record.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.

S.No.	Protocol Number	Product Name	Generic Name	Specific	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	SS Aspirin	2022 16:51:02	<a href="#">View</a>	<a href="#">View</a>
2	PS0023	Aspirin	Aspirin and codeine	SS Aspirin	2022 10:43:18	<a href="#">View</a>	<a href="#">View</a>

Go to page: 1 Row count: 10

**Comments \***: Approved


[Approve](#)
[Revert](#)

**Success**

Submitted successfully.

[OK](#)

LIMS 2.0
 All © rights Reserved PharmaSoftSol

	<b>Test Case Name</b>	Approve Schedule Test
	<b>Test Case Number</b>	LIMS-SM-008
	<b>Actual Result</b>	Pass


9. Click on User Name

SSNo : PSS-LIMS-SS-009

**PHARMA SOFT SOL**
**STABILITY MANAGEMENT**
LIMS APPROVER

### Schedule Approval

**Activity :**  
 . Select The Record.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.



Logged in as: **Protocol Modify Details**

S.No.	Protocol Number	Product Name	Generic Name	Specification Name	Customer	Registered By/On	Schedule Tests	Protocol	Protocol Modify Details
1	PS0010	Aspirin	Aspirin and codeine	SS Aspirin Spec	Dr Reddys	LIMS USER 27/06/2022 16:51:02	<a href="#">View</a>	<a href="#">View</a>	<a href="#">View</a>

<< < 1 > >> Go to page: 1 Row count: 10 Showing 1-1 of 1

**Comments \***

Maximum allowed characters should be 1500 only

[Approve](#)
[Revert](#)

**LIMS 2.0**
All © rights Reserved PharmaSoftSol

10. Click on Logout

SSNo : PSS-LIMS-SS-010

# LIMS



User Name

Password

Module

Select

[Login](#)
[Reset](#)