	<b>Test Case Name</b>	Reject Sample Acknowledgement
	<b>Test Case Number</b>	PSS-LIMS-025
	<b>Actual Result</b>	Pass

1. Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001





The login form is titled "LIMS" and features the PharmaSoftSol logo. It contains three input fields: "User Name" with the value "LIMS\_APPROVER", "Password" with masked characters "\*\*\*\*\*", and "Module" with a dropdown menu set to "STABILITY MANAGEMENT". Below the fields are "Login" and "Reset" buttons.

2. Click on Sample Acknowledgement

SSNo : PSS-LIMS-SS-002

**PHARMA SOFT**
**STABILITY MANAGEMENT**



LIMS APPROVER


### Sample Acknowledgment Approval

**Activity :**  
 . Select The Record and Click On view.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.

S.No	Protocol Number	Product Name	Specification Number	AR Number	Storage Condition	Stability Condition	Packing type	Batch Number	Start Date	Done By/On	Action
1	PS0014	Aspirin	SNSS 001	ASS220038	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	Product Batch No	21/06/2022	LIMS USER /21/06/2022 14:35:27	<a href="#">View</a>
2	PS0014	Aspirin	SNSS 001	ASS220041	30°C ± 2°C / 65%RH ± 5%RH	INTERMEDIATE	Box	Product Batch No	21/06/2022	LIMS USER /21/06/2022 14:40:00	<a href="#">View</a>
3	PS0024	Aspirin	SNSS 001	ASS220070	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	44	15/07/2022	LIMS USER /15/07/2022 11:12:29	<a href="#">View</a>

Go to page: 1 Row count: 10
Showing 1-3 of 3

LIMS 2.0
All © rights Reserved PharmaSoftSol

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Sample Acknowledgement
	<b>Test Case Number</b>	PSS-LIMS-025
	<b>Actual Result</b>	Pass

### 3. Select a Record And Click on View

SSNo : PSS-LIMS-SS-003

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS APPROVER

**Sample Acknowledgement**  
Activity :  
. Select The  
. Enter Com  
. Enter E-Sig  
**List of Samples**  

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0024  
Specification Number: SNSS 001  
Batch Number: 44  
Storage Conditions: 40°C ± 2°C / 75%RH ± 5%RH  
Qty Withdraw: 2 UOM: Gm  
Start Date: 15/07/2022  
Used For: Sampling  
Comments\*: Maximum allowed characters should be 1500  
Approve Revert

Product Name: Aspirin  
AR Number: ASS220070  
Packing Type: Box  
Interval: MONTH  
Total Qty Remained: 0 UOM: Gm  
View  
View  
View  
Showing 1-3 of 3

LIMS 2.0

### 4. Enter Comments

SSNo : PSS-LIMS-SS-004

**PHARMA SOFT**
**STABILITY MANAGEMENT**
LIMS APPROVER


**Sample Acknowledgement**  
Activity :  
. Select The  
. Enter Com  
. Enter E-Sig  
**List of Samples**  

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0024  
Specification Number: SNSS 001  
Batch Number: 44  
Storage Conditions: 40°C ± 2°C / 75%RH ± 5%RH  
Qty Withdraw: 2 UOM: Gm  
Start Date: 15/07/2022  
Used For: Sampling  
Comments\*: Rejected  
Approve Revert


Product Name: Aspirin  
AR Number: ASS220070  
Packing Type: Box  
Interval: MONTH  
Total Qty Remained: 0 UOM: Gm  
View  
View  
View  
Showing 1-3 of 3

LIMS 2.0

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Sample Acknowledgement
	<b>Test Case Number</b>	PSS-LIMS-025
	<b>Actual Result</b>	Pass

5. Click on Reject

SSNo : PSS-LIMS-SS-005



**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS APPROVER

**Sample Acknowledgement**

Activity :  
 . Select The  
 . Enter Com  
 . Enter E-Sig

**List of Samples**

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0024 Product Name: Aspirin

Specification Number: SNSS 001 AR Number: ASS220070

Batch Number: 44 Packing Type: Box

Storage Conditions: 40°C ± 2°C / 7

Qty Withdraw: 2 UOM: Gm

Start Date: 15/07/2022

Used For: Sampling

Comments\*: Rejected

**Validate E-Signature**

E-Signature Password\* :


Submit Cancel

Approve Revert

LIMS 2.0

6. Enter E-Signature Password

SSNo : PSS-LIMS-SS-006



**PHARMA SOFT** **STABILITY MANAGEMENT** LIMS APPROVER

**Sample Acknowledgement**

Activity :  
 . Select The  
 . Enter Com  
 . Enter E-Sig

**List of Samples**

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0024 Product Name: Aspirin

Specification Number: SNSS 001 AR Number: ASS220070

Batch Number: 44 Packing Type: Box

Storage Conditions: 40°C ± 2°C / 7

Qty Withdraw: 2 UOM: Gm

Start Date: 15/07/2022

Used For: Sampling

Comments\*: Rejected


**Validate E-Signature**

E-Signature Password\* :

Submit Cancel

Approve Revert

LIMS 2.0

 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Reject Sample Acknowledgement
	<b>Test Case Number</b>	PSS-LIMS-025
	<b>Actual Result</b>	Pass

7. Click on Submit

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT** **STABILITY MANAGEMENT** 88 LIMS APPROVER

Sample Acknowledgement

Activity :  
 . Select The  
 . Enter Com  
 . Enter E-Sig

List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0024

Product Name: Aspirin

Specification Number: SNSS 001

AR Number: ASS220070

Batch Number: 44

Packing Type: Box

Storage Conditions: 40°C ± 2°C / 75% RH ± 5%

Qty Withdraw: 2

UOM: Gm

Start Date: 15/07/2022

Used For: Sampling

Comments\*: Rejected

Approve Revert

**Validate E-Signature**

E-Signature Password\* : \*\*\*\*\*

Saving...

Submit Cancel

8. Click on OK button

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT** **STABILITY MANAGEMENT** 88 LIMS APPROVER

Sample Acknowledgement

Activity :  
 . Select The  
 . Enter Com  
 . Enter E-Sig

List of Samples

S.No	Protocol Number
1	PS001
2	PS001
3	PS002

Protocol Number: PS0024

Product Name: Aspirin

Specification Number: SNSS 001

AR Number: ASS220070

Batch Number: 44

Packing Type: Box

Storage Conditions: 40°C ± 2°C / 75% RH ± 5%

Qty Withdraw: 2

UOM: Gm

Start Date: 15/07/2022

Used For: Sampling


Comments\*: Rejected

Approve Revert

**Success**

Submitted successfully.


OK

	<b>Test Case Name</b>	Reject Sample Acknowledgement
	<b>Test Case Number</b>	PSS-LIMS-025
	<b>Actual Result</b>	Pass

9. Click on User Name


SSNo : PSS-LIMS-SS-009

**PHARMA SOFT SOL**
**STABILITY MANAGEMENT**


**LIMS APPROVER**

### Sample Acknowledgment Approval

**Activity :**  
 . Select The Record and Click On view.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.



S.No	Protocol Number	Product Name	Specification Number	AR Number	Storage Condition	Stability Condition	Packing type	Batch Number	Start Date	Done By/On	Logon	Action
1	PS0014	Aspirin	SNSS 001	ASS220038	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	Product Batch No	21/06/2022	LIMS USER / 21/06/2022 14:35:27		<a href="#">View</a>
2	PS0014	Aspirin	SNSS 001	ASS220041	30°C ± 2°C / 65%RH ± 5%RH	INTERMEDIATE	Box	Product Batch No	21/06/2022	LIMS USER / 21/06/2022 14:40:00		<a href="#">View</a>

Go to page:  Row count:  Showing 1-2 of 2

**LIMS 2.0**
All © rights Reserved PharmaSoftSol

10. Click on Logout

SSNo : PSS-LIMS-SS-010



## LIMS

User Name

Password

Module