	<b>Test Case Name</b>	Approve Protocol
	<b>Test Case Number</b>	LIMS-SM-005
	<b>Actual Result</b>	Pass

1.Enter the username, password, Select Module and click on login button


SSNo : PSS-LIMS-SS-001



The screenshot shows the LIMS login interface. At the top center, the word "LIMS" is displayed in yellow. Below it is the PharmaSoftSol logo. The login form contains three input fields: "User Name" with the value "LIMS\_APPROVER", "Password" with masked characters "\*\*\*\*\*", and "Module" with a dropdown menu set to "STABILITY MANAGEMENT". At the bottom of the form are two buttons: "Login" and "Reset".

2.Click on Protocol

SSNo : PSS-LIMS-SS-002



The screenshot shows the "Protocol Approval" screen in the LIMS application. The header includes the PharmaSoftSol logo, the text "STABILITY MANAGEMENT", and a user profile icon labeled "LIMS Approver". The main content area is titled "Protocol Approval" and contains instructions: "Activity : . Select The Record and Click On View & Approval. . Enter Comments and Click on Approve / Revert. . Enter E-Signature Password and Click on Submit." Below the instructions is a table titled "List of Protocols for Approval".


S.No.	Protocol Number	Product Name	Specification Name	Generic Name	Registered By/On	Request From	Action
1	PC0015	PARA65	SPPARACHEMICAL	Paracetamol 650 mg	LIMS User 31/08/2022 20:35:09	REVERT-RESUBMIT	<a href="#">View &amp; Approval</a>

Below the table, there are navigation controls: "<< < 1 > >>" and "Go to page 1" (with a dropdown), "Row count: 10" (with a dropdown), and "Showing 1-1 of 1".

The footer of the page displays "LIMS 2.0" on the left and "All © rights Reserved PharmaSoftSol" on the right.





 <b>PHARMA SOFT SOL</b> <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	<b>Test Case Name</b>	Approve Protocol
	<b>Test Case Number</b>	LIMS-SM-005
	<b>Actual Result</b>	Pass

7. Click on Submit

SSNo : PSS-LIMS-SS-007

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
👤 LIMS Approver ▾

Packing type	No	Batch size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Type
Blister	44	100	08/31/2022	08/31/2028	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH

Comments\*:

Approved

Validate E-Signature

E-Signature Password\*:

Saving...

Submit

Cancel

Approve

Revert

8. Click on OK button

SSNo : PSS-LIMS-SS-008

**PHARMA SOFT**
**STABILITY MANAGEMENT**
☰
👤 LIMS Approver ▾

Packing Type	Product Batch No	Batch Size	Manufactured Date	Expiry Date	Storage Condition	Orientation	Stability Condition	Interval Type	Int
Blister	44	100	08/31/2022	08/31/2028	40°C ± 2°C / 75%RH ± 5%RH	UP RIGHT	ACCELERATED	MONTH	

Comments\*:

Approved


Success

Submitted successfully.

OK

Approve


Revert

	<b>Test Case Name</b>	Approve Protocol
	<b>Test Case Number</b>	LIMS-SM-005
	<b>Actual Result</b>	Pass

9.Click on User Name


SSNo : PSS-LIMS-SS-009

**PHARMA SOFT SOL**
**STABILITY MANAGEMENT**

 LIMS Approver

### Protocol Approval

**Activity :**  
 . Select The Record and Click On View & Approval.  
 . Enter Comments and Click on Approve / Revert.  
 . Enter E-Signature Password and Click on Submit.



#### List of Protocols for Approval

S.No.	Protocol Number	Product Name	Specification Name	Generic Name	Registered By/On	Request From	Action
No data available!							


Row count: 10

LIMS 2.0
 All © rights Reserved PharmaSoftSol

10.Click on Logout

SSNo : PSS-LIMS-SS-010

# LIMS



User Name

Password

Module

Select

Login

Reset