
	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

1.Enter the username, password, Select Module and click on login button

SSNo : PSS-LIMS-SS-001



The screenshot shows the LIMS login interface. At the top, the 'LIMS' title is displayed in yellow. Below it, the PharmaSoftSol logo is visible. The login form contains three input fields: 'User Name' with the value 'LIMS_USER', 'Password' with masked characters '*****', and 'Module' with a dropdown menu set to 'STABILITY MANAGEMENT'. At the bottom of the form are two buttons: 'Login' and 'Reset'.

2.Click on Sample Acknowledgement


SSNo : PSS-LIMS-SS-002



The screenshot displays the 'Sample Acknowledgement' screen within the 'STABILITY MANAGEMENT' module. The header shows 'PHARMA SOFT' and 'STABILITY MANAGEMENT'. The user is logged in as 'LIMS USER'. The main content area has a tabbed interface with 'Sample Acknowledgement Details' selected. Below the tabs, there are instructions: '. Select Sample Acknowledgement, Click On Search, Select The Record and Click on Next.' and '. If Reverted Sample Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.' Two radio buttons are present: 'Sample Acknowledgement' (unselected) and 'Reverted Sample Acknowledgments by Approver' (selected). Below the radio buttons is a table titled 'List of Reverted Sample Acknowledgements by Approver'.

S.No	Protocol Id	Product Name	Specification Number	AR Number	Protocol Number	Storage Conditions	Stability Conditions	Packing Type	Batch Number	Start Date	RevertedBy/On	Comments
1	7	Aspirin	SNSS 001	ASS220013	PS0003	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	Product Batch No 1	11/06/2022	LIMS APPROVER 13/07/2022	view

Below the table, there is a pagination control showing 'Go to page: 1' and 'Row count: 10'. The status 'Showing 1-1 of 1' is displayed at the bottom right of the table area.

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

3. Click on Sample Acknowledgment

SSNo : PSS-LIMS-SS-003

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Sample Acknowledgment

Action Sample Acknowledgment Details

. Select Sample Acknowledgment, Click On Search, Select The Record and Click on Next.
. If Reverted Sample Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.

☒ Sample Acknowledgment
☐ Reverted Sample Acknowledgments by Approver

From Date : To Date:
Product : **Browse** Specification:

Search **Clear**

Samples Loaded														
S.No	Product	Specification Number	AR Number	Protocol Number	Storage Condition	Stability Condition	Orientation	Packing Type	Batch Number	Start Date	Total Qty Remained	Qty Remained UOM	Interval	Is deviated
No data available!														

Go to page: Row count:

IMS 2.0
All © rights Reserved PharmaSoftSol

4. Select a Record

SSNo : PSS-LIMS-SS-004

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Sample Acknowledgment

Action Sample Acknowledgment Details

. Select Sample Acknowledgment, Click On Search, Select The Record and Click on Next.
. If Reverted Sample Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.


☒ Sample Acknowledgment
☐ Reverted Sample Acknowledgments by Approver

From Date : To Date:
Product : **Browse** Specification:

Search **Clear**

Samples Loaded														
S.No	Product	Specification Number	AR Number	Protocol Number	Storage Condition	Stability Condition	Orientation	Packing Type	Batch Number	Start Date	Total Qty Remained	Qty Remained UOM	Interval	Is deviated
1	Aspirin	SNSS 001	A5S220014	PS0003	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	UP RIGHT	Box	Product Batch No 1	11/06/2022	43	Mg	3 MONTH	No
					40°C ± 2°C /				Product					

IMS 2.0
All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

5. Click on Next

SSNo : PSS-LIMS-SS-005

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Sample Acknowledgment

[Action](#)
[Sample Acknowledgment Details](#)

Click On Search, Select The Record and Click on Next.
Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.

☐ Acknowledgment
 ☐ Reverted Sample Acknowledgments by Approver

To Date:

Specification:

ion	AR Number	Protocol Number	Storage Condition	Stability Condition	Orientation	Packing Type	Batch Number	Start Date	Total Qty Remained	Qty Remained UOM	Interval	Is deviated
	ASS220014	PS0003	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	UP RIGHT	Box	Product Batch No 1	11/06/2022	43	Mg	3 MONTH	No
			40°C ± 2°C /				Product					

IMS 2.0
All © rights Reserved PharmaSoftSol

6. Enter Quantity

SSNo : PSS-LIMS-SS-006

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Sample Acknowledgment

[Action](#)
[Sample Acknowledgment Details](#)

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :

UOM :


Approval From:

Used For:

Maximum allowed characters should be 100

Select an Action*:
☐ Acknowledge
☐ Terminate

IMS 2.0
All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

7.Select UOM
SSNo : PSS-LIMS-SS-007

PHARMA SOFT

STABILITY MANAGEMENT

LIMS USER

Sample Acknowledgment

Action

Sample Acknowledgment Details

• Enter Quantity withdraw and select UOM from drop down.click on Select to select the Approver.

• If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish

• Enter E- Signature Password and Click On Finish.

Quantity Withdraw :

2

UOM :

Gm

Approval From:

Select

Used For:

Maximum allowed characters should be 100

Select an Action*:

☐ Acknowledge

☐ Terminate

I IMS 2.0

All © rights Reserved PharmaSoftSol

8.Enter Used For
SSNo : PSS-LIMS-SS-008

PHARMA SOFT

STABILITY MANAGEMENT

LIMS USER

Sample Acknowledgment

Action

Sample Acknowledgment Details

• Enter Quantity withdraw and select UOM from drop down.click on Select to select the Approver.

• If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish

• Enter E- Signature Password and Click On Finish.

Quantity Withdraw :

2

UOM :

Gm

Approval From:

Select

Used For:

Sampling


Select an Action*:

☐ Acknowledge

☐ Terminate


I IMS 2.0

All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

11.Select a Record

SSNo : PSS-LIMS-SS-011




The screenshot shows the 'STABILITY MANAGEMENT' application interface. On the left, there is a sidebar with 'Sample Acknowledgement' and 'Sample Acknowledgement Details' tabs. The main area is titled 'Locations' and contains a tree view with 'organization' and 'Plant-1' (which has sub-items 'QA', 'QC', and 'PF'). To the right of the tree view are input fields for 'First Name' and 'Last Name', and 'Search' and 'Clear' buttons. Below this is a 'Users' table with columns 'First Name' and 'Last Name'. The table contains three rows: 'LIMS USER', 'LIMS REVIEWER', and 'LIMS APPROVER'. The 'LIMS APPROVER' row is highlighted. At the bottom of the table are 'Select' and 'Cancel' buttons. The footer shows 'LIMS 2.0' and 'All © rights Reserved PharmaSoftSol'.

12.Click on Select

SSNo : PSS-LIMS-SS-012



The screenshot shows the 'STABILITY MANAGEMENT' application interface. The main area is titled 'Sample Acknowledgement' and contains a 'Sample Acknowledgement Details' tab. Below the tab are instructions: 'Enter Quantity withdraw and select UOM from drop down.click on Select to select the Approver.', 'If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish', and 'Enter E- Signature Password and Click On Finish.'. The form has input fields for 'Quantity Withdraw' (value 2), 'UOM' (value Gm), and 'Approval From' (value LIMS APPROVER). There is a 'Select' button next to 'Approval From'. Below these is a 'Used For' text area with the value 'Sampling'. At the bottom, there is a 'Select an Action' section with radio buttons for 'Acknowledge' and 'Terminate'. The footer shows 'LIMS 2.0' and 'All © rights Reserved PharmaSoftSol'.

 PHARMA SOFT SOL <small>COMMITTED TO PAPERLESS COMPLIANCE</small>	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

13. Click on Acknowledge

SSNo : PSS-LIMS-SS-013

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Sample Acknowledgment

Action
Sample Acknowledgment Details

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :
UOM :
Approval From:

Used For:

Select an Action*: ☒ Acknowledge ☐ Terminate

LIMS 2.0
All © rights Reserved PharmaSoftSol

14. Click on Finish

SSNo : PSS-LIMS-SS-014

PHARMA SOFT
STABILITY MANAGEMENT
LIMS USER

Sample Acknowledgment

Action
Sample Acknowledgment Details

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.


Quantity Withdraw :
UOM :
Approval From:

Used For:

Select an Action*: ☒ Acknowledge ☐ Terminate

Validate E-Signature
E-Signature Password*:

LIMS 2.0
All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

15. Enter E-Signature Password

SSNo : PSS-LIMS-SS-015

PHARMA SOFT **STABILITY MANAGEMENT** LIMS USER

Sample Acknowledgment

Action **Sample Acknowledgment Details**

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :

Used For:

Select an Action*: ☒ Acknowledgment

Validate E-Signature

E-Signature Password*:

LIMS APPROVER

LIMS 2.0 All © rights Reserved PharmaSoftSol

16. Click on Submit

SSNo : PSS-LIMS-SS-016

PHARMA SOFT **STABILITY MANAGEMENT** LIMS USER

Sample Acknowledgment

Action **Sample Acknowledgment Details**

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw :

Used For:


Select an Action*: ☒ Acknowledgment

Validate E-Signature

E-Signature Password*:

LIMS APPROVER

LIMS 2.0 All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

17. Click on OK button

SSNo : PSS-LIMS-SS-017

PHARMA SOFT **STABILITY MANAGEMENT** LIMS USER

Sample Acknowledgment

Action **Sample Acknowledgment Details**

- Enter Quantity withdraw and select UOM from drop down. click on Select to select the Approver.
- If select Acknowledgment then enter Used for comments , or if select Terminate, then enter Comment if Terminate, click finish
- Enter E- Signature Password and Click On Finish.

Quantity Withdraw : UOM :

Used For:

Select an Action*: ☒ Acknowledge ☐ Terminate

Success

Submitted successfully.

LIMS 2.0 All © rights Reserved PharmaSoftSol

18. Click on User Name

SSNo : PSS-LIMS-SS-018

PHARMA SOFT **STABILITY MANAGEMENT** LIMS USER

Sample Acknowledgment

Action **Sample Acknowledgment Details**


- . Select Sample Acknowledgment, Click On Search, Select The Record and Click on Next.
- . If Reverted Sample Acknowledgments by Approver is Required, Then Select Reverted Sample Acknowledgments by Approver.

☐ Sample Acknowledgment ☒ Reverted Sample Acknowledgments by Approver

S.No	Protocol Id	Product Name	Specification Number	AR Number	Protocol Number	Storage Conditions	Stability Conditions	Packing Type	Batch Number	Start Date	RevertedBy/On	Comments
1	7	Aspirin	SNSS 001	ASS220013	PS0003	40°C ± 2°C / 75%RH ± 5%RH	ACCELERATED	Box	Product Batch No 1	11/06/2022	LIMS APPROVER 13/07/2022	<input type="button" value="view"/>

<< < 1 > >> Go to page: Row count: Showing 1-1 of 1

LIMS 2.0 All © rights Reserved PharmaSoftSol

	Test Case Name	Create Sample Acknowledgement
	Test Case Number	LIMS-SM-013
	Actual Result	Pass

19. Click on Logout

SSNo : PSS-LIMS-SS-019

